



**SPECTRUM ANALYTICAL, INC.**  
 Featuring  
**HANIBAL TECHNOLOGY**

**Sample Container Request**

Date Ordered: \_\_\_/\_\_\_/\_\_\_ Company Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_ Zip: \_\_\_\_\_

Date Needed: \_\_\_/\_\_\_/\_\_\_ Contact Person: \_\_\_\_\_

STOCK	CASE QTY	# OF INDIVIDUAL NEEDED (DO NOT ORDER BY THE CASE)	*PRESERVATIVE
<b>Containers for Aqueous Samples:</b>			
40ml Amber Vials (VOC)	72		HCl
40ml Amber Vials (EDB by 504)	72		Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>
40ml Amber Vials (TOC)	72		H <sub>3</sub> PO <sub>4</sub>
1L Amber Glass	12		
½ Gallon Plastic	6		
1L Plastic	12		
500ml Plastic	24		
250ml Plastic	24		
100ml Sterile Plastic (Coliforms)	1		Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>
<b>Containers for Solid Samples:</b>			
8oz. Amber Soil Jar	12		
4oz. Amber Soil Jar	24		
40ml Clear Vial	72		Deionized H <sub>2</sub> O
40ml Clear Vial	72		CH <sub>3</sub> OH
40ml Clear Vial	72		Not preserved
<b>Containers for Air Samples:</b>			
1L Tedlar Bag	10		
Summa Canisters/Regulators		<i>Please use Air Sampling Media Form</i>	
<b>Other Supplies:</b>			
Chains of Custody			
Sample Labels			
Trip Blanks			
Coolers			S – M – L

\* Please note laboratory preservation if needed.

Received by: \_\_\_\_\_ Date: \_\_\_\_\_