



SPECTRUM ANALYTICAL, INC.  
*Featuring*  
HANIBAL TECHNOLOGY

# SAMPLING PROTOCOL FOR ACIDITY & ALKALINITY

## COLLECTION OF AQUEOUS SAMPLES FOR ACIDITY AND ALKALINITY

For more information, please feel free to contact our Quality Service Department at (800) 789-9115.

Standard Methods for the Examination of Water and Wastewater sampling and storage protocols advise that waste samples may be subject to microbial action and loss or gain of CO<sub>2</sub> or other gases when exposed to air. To help maintain sample integrity, sample containers for acidity and/or alkalinity should be filled completely and capped tightly immediately upon collection.